


## EXPERIMENT: ANALYSIS OF A FREELY FALLING BODY

### Part I: dependence of speed and position on time.

Primary OBJECTIVE: Investigate whether the velocity of the falling body increases linearly with time and calculate a value for  $g$ , the acceleration due to gravity. Verify how the distance varies as a function of time

### APPARATUS:

Behr Freefall System ,Spark Timing System



### THEORY

In this experiment, a brass cylinder is dropped and a record of its fall is made. Before the measurement, the cylinder is suspended at the top of the stand with the help of an electromagnet. When the electromagnet is turned off, the cylinder starts to fall. Simultaneously, the spark timer starts to send short high-voltage pulses on two wires, stretched along the cylinder's path.

At the time of each pulse a spark goes through the wires and the cylinder, leaving a mark on the special paper tape. The time interval between two adjacent sparks  $\tau$  is  $1/60$  of a second.

Measuring the distances between any two marks  $\Delta y$  and knowing the times between the corresponding sparks  $\Delta t$ , it is possible to calculate the average velocity of the cylinder in this interval, using the formula

$$v = \frac{\Delta y}{\Delta t}. \quad (1)$$

If  $\Delta t$  is small enough, we can assume that the velocity at any instant within this interval is approximately equal to this average velocity.

For the motion of the body with the constant acceleration  $a$ , the velocity  $v$  is given by

$$v = at + v_o, \quad (2)$$

where  $v_o$  is the velocity of the cylinder at  $t = 0$ . Since the brass cylinder is falling freely,

$$a = -g, \quad (3)$$

where  $g=9.81 \text{ m/sec}^2$  is the magnitude of the gravity acceleration, it follows from (2) that

$$v = v_o - gt \quad (4)$$

Thus  $g$  can be determined from a plot of  $v$  vs.  $t$ , and compared with the known value. The position of the cylinder, denoted by  $y$ , as a function of time, is given by the standard equation for an object which is undergoing constant acceleration. If it starts to fall from an initial height  $y_o$  with an initial velocity  $v_o$ , then this equation looks like

$$y = y_o + v_o t - \frac{1}{2} g t^2 \quad (5)$$

## PROCEDURE

Using the free fall apparatus, drop the brass cylinder and record on spark tape the location of the falling cylinder after a series of equally spaced time intervals,  $\tau = 1/60 \text{ sec}$ .

The quantities needed to analyze the motion are the position ( $y$ ), velocity ( $v$ ) and time ( $t$ ) of the points on your spark tape. The choices of time,  $t = 0$  and height  $y = 0$  are arbitrary and do not necessarily refer to your first or last data points. However, as your object falls the position  $y$  must decrease (becoming negative if necessary).

After you have performed the experiment, tape the paper strip to your lab table. Label the points on your tape, starting with the last point as #31 and label them in descending order. Make sure that point #31 is the last data point and not the point where the brass cylinder is stuck in the putty and simply sparking. Point #1 will not necessarily correspond to the very first point, but the very first points are somewhat ambiguous anyway. Put a ruler on your tape such that  $y=0.0 \text{ cm}$  corresponds to point #31 and  $y$  increases as you move towards point #1. Measure the locations of each marker with this ruler and write the  $y$  position on the tape next to each marker.

After completing this, transfer the  $y$  positions into the spreadsheet. Calculate the velocity at all points, except the first and the last, during the fall using the spreadsheet. Be sure you do at least one sample calculation. The instantaneous velocities at these points are

approximated by finding the average velocity at each of your measurements  $v = \frac{\Delta y}{\Delta t}$ ,

where  $\Delta y$  is the distance between the preceding spark and the following spark and  $\Delta t$  is the time interval between these sparks. More details are given in the Appendix "Formulas, definitions and errors for the Free Fall experiment".

Transfer the columns of time,  $y$  positions and velocities from your spreadsheet to Kaleidagraph for graphing.

***Prepare the following graphs:***

I -- a graph showing  $v$  vs.  $t$  with the slope and intercept clearly indicated and calculated. From the slope, determine the gravitational acceleration  $-g$  in  $\text{cm}/\text{sec}^2$ .

II-- a graph showing the  $y$  position versus time.

*Questions to be answered:*

- 1) When finding the velocities why should you use sparks on either side of the point at which you are calculating the speed by  $v = \frac{\Delta y}{\Delta t}$ ?
- 2) What is the intercept (value of velocity at  $t=0$  sec) determined from graph I and what does it mean?
- 3) From graph I find the time at which  $v=0$   $\text{cm}/\text{sec}$ .
- 4) Determine the slope from graph I and from that the gravitational acceleration  $-g$  in  $\text{cm}/\text{sec}^2$ .
- 5) What equation describes graph II?
- 6) When the initial velocity is zero, what would you plot to make graph II linear:  $y^2$  vs.  $t$  or  $y^2$  vs.  $t^2$  or  $y$  vs.  $t^2$ ?

**CHECKLIST:**

Your lab report should include the following:

- 1) A spreadsheet with point numbers and  $\Delta y$ 's recorded + sample calculations.
- 2) Plot of the height vs. time
- 3) Plot of velocity vs. time with slope (in  $\text{cm}/\text{sec}^2$ )
- 4) Interpretation of the two plots
- 5) Answers to questions
- 6) One member of each group should turn in a spark tape record of the freefall
- 7) Since You will use the obtained results in the next lab, save the spreadsheet in the folder "Results/Your section". Do not forget to take the manual for the next lab!

## Formulas, definitions and errors for the Free Fall Experiment

### Definitions:

In this experiment you measure the position of a falling mass  $-m-$  at fixed time intervals. The fixed time interval is determined by a high voltage spark source. Read off the time between sparks ( $-\tau-$ ) from the setting on the spark source.

You will measure the positions at each spark as  $y_1, y_2, y_3, y_4,$  etc. in centimeters [cm].

These positions will be referred to as  $y_i$ .

To measure the speed at point  $-i-$  (for example  $i=6$ ) you determine:

$$\Delta y_i = y_{i+1} - y_{i-1} \text{ in our particular example: } \Delta y_6 = y_7 - y_5$$

from your measurements. On your data sheet this is labeled as  $\Delta y(i)$ . The speed at point  $-i-$  is now determined by dividing the distance  $\Delta y_i$  ( $=\Delta y(i)$ ) by the time elapsed between the two points:  $2\tau$ . So the speed  $V_{y_i}$  is given by:

$$V_{y_i} = \frac{\Delta y_i}{\Delta t} = \frac{\Delta y_i}{2\tau} \text{ [cm/sec]}$$

### The errors for these quantities are defined as follows:

For each measured  $y_i$  you assign an error based on how accurate you can measure that point. This error is called  $\delta y$ . This error determines all other errors in this lab. For the following formulas it is assumed that the error in  $-\tau-$  and  $-m-$  are zero. Here are the definitions:

The error in  $\Delta y$  at each point  $-i-$  is the same and is given by  $\delta(\Delta y) = 2\delta y$

The error in the speed at each point  $-i-$  is:

$$\delta(V_y) = V_y \frac{\delta(\Delta y)}{\Delta y} = V_y \frac{2\delta y}{\Delta y} = \frac{V_y}{\Delta y} 2\delta y = \frac{2\delta y}{2\tau} = \frac{\delta y}{\tau}$$

Free Fall 1\_empty

Free Fall with Behr apparatus				
Partner:			Name:	
			Section:	
			value	units
Falling mass m=			0.01	[grams]
Error in mass is assumed to be zero.				
Time Interval between sparks:			0.0167	[seconds]
Gravitational acceleration -g-:			980	[cm/sec**2]
Error in -y- : (y)			xxx	[cm]
Error in - y- : ( y)= 2 (y)			xxx	[cm]
Error in -Vy- : (Vy)= ( y)*Vy/ y			##	[cm/sec]
			Recall:	
			Blue = measured values	
			Yellow = calculated values	
$y(i)=y(i+1) - y(i-1)$				
Dot	Time	y(i)	y(i)	Vy(i)
i	[sec]	[cm]	[cm]	[cm/sec]
1	0.000			
2	0.017			
3	0.033			
4	0.050			
5	0.067			
6	0.083			
7	0.100			
8	0.117			
9	0.133			
10	0.150			
11	0.167			
12	0.183			
13	0.200			
14	0.217			
15	0.233			
16	0.250			
17	0.267			
18	0.283			
19	0.300			
20	0.317			
21	0.333			
22	0.350			
23	0.367			
24	0.383			
25	0.400			
26	0.417			
27	0.433			
28	0.450			
29	0.467			
30	0.483			
31	0.500			